

OXYGEN MASK HIGH CONCENTRATION



MAIN SPECIFICATIONS

- Mask for high throughput oxygen administration.
- Transparent for better viewing.
- Single use.
- Metal clip to better fit nose.
- Bag reservoir for better control of the administration.
- Lateral filters.
- Elastic for fixing the mask to the head.
- Tube 2 m.

Code	Description	Selling unit
39011048	ADULT OXYGEN MASK HIGH CONCENTRATION WITH TUBING AND BAG	CX.100
39011049	PEDIATRIC OXYGEN MASK HIGH CONCENTRATION WITH TUBING AND BAG	CX.100