

ASKINA BIOFILM TRANSPARENT



MAIN SPECIFICATIONS

- Askina Biofilm transparent is a thin sterile transparent hydrocolloid dressing for moist wound treatment with the capacity to absorb residual exudate.
- Composition: Sodium-carboxymethylcellulose, Polyisobutylene and Polyurethane film.
- Thin design that fits any body contours.
- Transparent, allows visual control of the wound.
- Easy to apply.
- Stays in place during showering.
- Management of low exuding superficial chronic wounds (pressure ulcers, leg ulcers), superficial acute wounds (1st and 2nd degree burns, abrasions) and post operative wounds.

Code	Description	Selling unit
01001003	ASKINA BIOFILM TRANSPARENT 10X10CM	CX. 10
F72092	ASKINA BIOFILM TRANSPARENT 20X20CM	CX. 5
F72091	ASKINA BIOFILM TRANSPARENT 15X15CM	CX. 5
F72095	ASKINA BIOFILM TRANSPARENT 5X20 CM	CX. 10