

# GRASSOLIND PARAFFIN COTTON DRESSING



## MAIN SPECIFICATIONS

- Grassolind is an open-weave cotton dressing, impregnated with a non-medicated paraffin with no sensitising effects.
- Grassolind is sufficiently porous to allow an easy passage of exudate from the wound surface, preventing maceration of surrounding tissue.
- The dressing is ideal for treatment of large surface wounds where low adherence is essential. The dressing can be kept on the wound for up to 7 days, subject to the amount of exudate. Pads without active substances of cottonknitted fabric batter. permeable to air and in the secretions, impregnated with a neutral ointment in a large amount in order not removed even in intense exudation. Readily cut to the desired dimensions.
- The ointment contains:
  - white petrolatum, ester of diglycerides, mono and dicarboxylic fatty acids, synthetic wax.
- uses
- For a general treatment of wounds, in particular for maintaining the elasticity of the surface and the edges of the wound, e.g. for extensive abrasions or cuts and burns or burns from scalding hot water, plastic and cosmetic surgery in exporting nail in phimosis surgery etc. To cover the donor and recipient region in skin transplants and to hold partial thickness skin grafts. They are also

indicated for dermatology and for patients with sensitive skin or sensitive to medications.

Code	Description	Selling unit
01046038	GRASSOLIND 7,5X10CM	CX. 10
01046039	GRASSOLIND 10X10 CM	CX. 10
01046044	GRASSOLIND 10X20 CM	CX. 30
4993102 (H)	GRASSOLIND 5X5 CM	CX. 10