

# THREE WAY STOPCOCK



## MAIN SPECIFICATIONS

- Rotating male luer lock and 2 female luer lock connectors.
- Full transparent / polycarbonate stopcock holder.
- 360° rotatable without limitation.
- Three way stopcock with extension tubing of various lengths.
- Available in Plain and With Extension Tubing, with or without injection valve.
- With protective cap and Luer-lock safety lock.

Code	Description	Selling unit
09011003	THREE WAY STOPCOCK	CX. 50
09011002	THREE WAY STOPCOCK WITH EXTENSION TUBING 10 CM	CX. 50
09085009	THREE WAY STOPCOCK WITH EXTENSION TUBING 25 CM LUER LOCK	CX. 50