

JELONET COTTON DRESSING



MAIN SPECIFICATIONS

- Composed of knitted cotton fabric impregnated with paraffin.
- Soothes and protects the wound.
- It allows free passage of exudate.
- It can be used in combination with a secondary dressing and antiseptics.
- Suitable for: burns, leg ulcers, wounds with loss of skin, lacerations and abrasions.
- They are provided sterile in a variety of sizes.
- Soothes and protects the wound and allows free passage of viscous exudate.
- It has a soothing low tack which allows the wound to drain freely into the absorbent secondary dressing.
- Retains shape without fraying, open woven gauze dressings made have connected wires to be cut not fray.
- It can be combined with topical medication of choice.
- It does not have medicine and is ideal for topical use with antiseptic or antibiotic of choice.
- Protects injury preventing adhesion of the dressing shield and preventing the dehydration of the granulation tissue.
- It is indicated for the treatment of skin lesions clean, protecting and preventing rupture of the granulation tissue on dressing change. Suitable for frequent changes.

| Code | Description | Selling unit |
|----------|------------------|--------------|
| 01003014 | JELONET 5X5 CM | CX. 50 |
| 01003013 | JELONET 10X10 CM | CX.10 |