

# COSMOPOR DRESSING INTRAVENOUS



## MAIN SPECIFICATIONS

- Self-adhesive I.V. dressing for the fixation of cannulas, non-woven fabric support material, with an absorbent pad with a special coating not adhering to the wound.
- Additional separate dressing for the protection of pressure points caused by the cannula.
- Good skin tolerance due to hypoallergenic synthetic adhesive.

Code	Description	Selling unit
01046021	COSMOPOR IV DRESSING POSTOPERATIVE 6X8 CM	CX. 50