

# DIAGNOSTIC GLOVES LATEX WITH POWDER PROHAND



## MAIN SPECIFICATIONS

- Diagnostic sterile latex gloves
- Glove type: diagnostic, pre-powdered, sterile
- Raw material: natural rubber latex
- Powdering agent: corn starch (meal)
- Quantity in an individual packaging: 1 pair, 50 pairs, 300 pairs
- Surface: smooth
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: S, M, L

Code	Description	Selling unit
-	DIAGNOSTIC GLOVES LATEX WITH POWDER PROHAND S	CX. 50
--	DIAGNOSTIC GLOVES LATEX WITH POWDER PROHAND M	CX. 50
---	DIAGNOSTIC GLOVES LATEX WITH POWDER PROHAND L	CX. 50



[www.factormed.pt](http://www.factormed.pt)