

EXAM GLOVES LATEX WITHOUT POWDER

DERMAGEL



MAIN SPECIFICATIONS

- Powder-free latex gloves, diagnostic and protective
- Glove type: diagnostic, protective, non-sterile
- Raw material: natural rubber latex
- Powdering agent: none - internal polymer layer
- Quantity in an individual packaging: 100 pieces
- Surface: textured
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES LATEX WITHOUT POWDER DERMAGEL XS	CX. 100
--	EXAM GLOVES LATEX WITHOUT POWDER DERMAGEL S	CX. 100
---	EXAM GLOVES LATEX WITHOUT POWDER DERMAGEL M	CX. 100

----	EXAM GLOVES LATEX WITHOUT POWDER DERMAGEL L	CX. 100
------	---	---------

-----	EXAM GLOVES LATEX WITHOUT POWDER DERMAGEL XL	CX. 100
-------	--	---------
