

EXAM GLOVES LATEX WITH POWDER

SANTEX



MAIN SPECIFICATIONS

- Pre-powdered latex gloves, diagnostic and protective
- glove type: diagnostic, protective, non-sterile
- raw material: natural rubber latex
- powdering agent: corn starch (meal)
- quantity in an individual packaging: 100 pieces
- surface: smooth/textured
- cuff finishing: evenly rolled edge
- shape: universal, fitting the right and left hand
- sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES LATEX WITH POWDER SANTEX XS	CX. 100
--	EXAM GLOVES LATEX WITH POWDER SANTEX S	CX. 100
---	EXAM GLOVES LATEX WITH POWDER SANTEX M	CX. 100

----	EXAM GLOVES LATEX WITH POWDER SANTEX L	CX. 100
-----	EXAM GLOVES LATEX WITH POWDER SANTEX XL	CX. 100