

EXAM GLOVES LATEX WITH POWDER

TEXIDENT MINT



MAIN SPECIFICATIONS

- Mint-scented pre-powdered latex gloves, diagnostic and protective, recommended for dentists
- Glove type: diagnostic, protective, non-sterile, mint-scented
- Raw material: natural rubber latex
- Powdering agent: corn starch (meal)
- Quantity in an individual packaging: 100 pieces
- Surface: micro-textured
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES LATEX WITH POWDER TEXIDENT MINT XS	CX. 100
--	EXAM GLOVES LATEX WITH POWDER TEXIDENT MINT S	CX. 100
---	EXAM GLOVES LATEX WITH POWDER TEXIDENT MINT M	CX. 100

----	EXAM GLOVES LATEX WITH POWDER TEXIDENT MINT L	CX. 100
------	---	---------

-----	EXAM GLOVES LATEX WITH POWDER TEXIDENT MINT XL	CX. 100
-------	--	---------
