

SYRINGE INJ/LIGHT (WITH NEEDLE, ECCENTRIC CONE)



MAIN SPECIFICATIONS

- PP barrell with graduated scale
- PP plunger
- Latex free gasket
- Eccentric Luer cone
- Triple bevelling stainless steel needle AISI 304
- PP needle cap
- Syringe type: standard syringes
- Safety mechanism: no

Code	Description	Selling unit
-	SYRINGE INJ/LIGHT (WITH NEEDLE, ECCENTRIC CONE) 10 ML NEEDLE 0,8X38 MM	CX. 100
--	SYRINGE INJ/LIGHT (WITH NEEDLE, ECCENTRIC CONE) 20 ML AGULHA 1,1X38 MM	CX. 50
---	SYRINGE INJ/LIGHT (WITH NEEDLE, ECCENTRIC CONE) 20 ML NEEDLE 0,9X38 MM	CX. 50

SYRINGE INJ/LIGHT (WITH NEEDLE, ECCENTRIC CONE) 20 ML
NEEDLE 0,8X38 MM

CX. 50
