


NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM



MAIN SPECIFICATIONS

- Stainless steel needle AISI 304 with triple bevelling, diameter and length according to ISO 9626
- PP needle protection cap
- Low memory 30 cm PVC tube or phthalates free (DBP,BBP,DEHP,DNOP, e DIDP)
- PP luer lock connector
- Type of needle: Safety scalp vein sets
- Safety mechanism: yes

Code	Description	Selling unit
-	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 18G - 1,2x19 MM PINK	CX. 100
--	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 19G - 1,1x19 MM PINK	CX. 100
---	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 20G - 0,9x19 MM YELLOW	CX. 100
----	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 21G - 0,8x19 MM GREEN	CX. 100

-----	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 22G - 0,7x19 MM BLACK	CX. 100
-----	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 23G - 0,6x19 MM BLUE	CX. 100
-----	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 25G - 0,5x19 MM ORANGE	CX. 100
-----	NEEDLE FLY/SAFE SCALP VEIN SET WITH SAFETY MECHANISM 27G - 0,4x19 MM GREY	CX. 100