

EXTENSIBLE FOLDING SCREEN SILENTIA



MAIN SPECIFICATIONS

- The screen can be mounted on a rotated carrier, thus allowing greater mobility.
- This rotated stand is suitable for extensible screens of various lengths: from 5 to 15 sections.
- Height: 1.45 m, 1.65 m or 1.85 m.
- Product available in several colours and patterns.
- Its patented design takes up minimal floor space while offering maximum stability.
- The extensible screens are compatible with each other and include the EasyClick function, so you can easily move a screen from a wall to a rotated stand - or combine a fixed screen and an extensible to take full advantage of all the possibilities the system offers .
- The screens have flat surfaces, resistant to disinfectants normally used in health care. The included wheels are easily removable and washable. The products meet the strict hygiene requirements, all of their surfaces being easily accessible for cleaning - this is EasyClean.
- The patented hinges have an integrated memory that allows for easy handling as well as a good prevention of cross-infections by having only a single point of contact: its handle with ergonomic design.
- Mobility with just one click. An extensible wall-mounted screen can be easily removable to a rotated bracket thus allowing its use with other components of the system. This is a complete system, of great versatility and of easy adaptation to the most diverse needs.
-

The wall mounting column of the extensible screen is usually on the left side, which implies that the screen is folded to the right of the fixing rail. For special situations, such as the proximity of a door, the screen can be ordered with the fixing system, the opposite side.

Code	Description	Selling unit
6105	FOLDING SCREEN OF 5 SECTIONS LENGTH 1.25M, HEIGHT 1.45M, WEIGHT 8KG	CX. 1
6205	FOLDING SCREEN OF 5 SECTIONS LENGTH 1.25M, HEIGHT 1.65M, WEIGHT 9.5KG	CX. 1
6305	FOLDING SCREEN OF 5 SECTIONS LENGTH 1.25M, HEIGHT 1.85M, WEIGHT 11.5KG	CX. 1
6107	FOLDING SCREEN OF 7 SECTIONS LENGTH 1.75M, HEIGHT 1.45M, WEIGHT 12KG	CX. 1
6207	FOLDING SCREEN OF 7 SECTIONS LENGTH 1.75M, HEIGHT 1.65M, WEIGHT 14KG	CX. 1
6307	FOLDING SCREEN OF 7 SECTIONS LENGTH 1.75M, HEIGHT 1.85M, WEIGHT 16KG	CX. 1
6109	FOLDING SCREEN OF 9 SECTIONS LENGTH 2.25M, HEIGHT 1.45M, WEIGHT 15.5KG	CX. 1
6209	FOLDING SCREEN OF 9 SECTIONS LENGTH 2.25M, HEIGHT 1.65M, WEIGHT 18KG	CX. 1
6309	FOLDING SCREEN OF 9 SECTIONS LENGTH 2.25M, HEIGHT 1.85M, WEIGHT 20.5KG	CX. 1
6111	FOLDING SCREEN OF 11 SECTIONS LENGTH 2.75M, HEIGHT 1.45M, WEIGHT 18.5KG	CX. 1
6211	FOLDING SCREEN OF 11 SECTIONS LENGTH 2.75M, HEIGHT 1.65M, WEIGHT 21.5KG	CX. 1
6311	FOLDING SCREEN OF 11 SECTIONS LENGTH 2.75M, HEIGHT 1.85M, WEIGHT 24.5KG	CX. 1
6213	FOLDING SCREEN OF 13 SECTIONS LENGTH 3.25M, HEIGHT 1.65M, WEIGHT 25.5KG	CX. 1
6113	FOLDING SCREEN OF 13 SECTIONS LENGTH 3.25M, HEIGHT 1.45M, WEIGHT 22KG	CX. 1
6313	FOLDING SCREEN OF 13 SECTIONS LENGTH 3.25M, HEIGHT 1.85M, WEIGHT 29KG	CX. 1
6115	FOLDING SCREEN OF 15 SECTIONS LENGTH 3.75M, HEIGHT 1.45M, WEIGHT 25.5KG	CX. 1
6215	FOLDING SCREEN OF 15 SECTIONS LENGTH 3.75M, HEIGHT 1.65M, WEIGHT 29.5KG	CX. 1

6315	FOLDING SCREEN OF 15 SECTIONS LENGTH 3.75M, HEIGHT 1.85M, WEIGHT 33.5KG	CX. 1
6514	WHEELED SUPPORT HEIGHT	CX. 1
6518	WHEELED SUPPORT HEIGHT 1.65M	CX. 1
6520	WHEELED SUPPORT HEIGHT 1.85M	CX. 1