

# SURGICAL STERILIZED GLOVE



## MAIN SPECIFICATIONS

- Natural latex surgical gloves with low content of protein, anatomically shaped, micro-roughened surface that prevents sliding of surgical instruments even in contact with liquids.
- Excellent sensitivity, double wall giving high protection for the physician and patient.
- In powdered gloves, the powder used is corn starch, perfectly inert material.
- Sterilized by gamma rays.

Code	Description	Selling unit
13011018	LUVA CIRÚRGICA ESTÉRIL COMFIT 6	CX.50
13011019	LUVA CIRÚRGICA ESTÉRIL COMFIT 6,5	CX.50
13011020	LUVA CIRÚRGICA ESTÉRIL COMFIT 7	CX.50
13011021	LUVA CIRÚRGICA ESTÉRIL COMFIT 7,5	CX.50
13011022	LUVA CIRÚRGICA ESTÉRIL COMFIT 8	CX.50
13011023	LUVA CIRÚRGICA ESTÉRIL COMFIT 8,5	CX.50