

# SYRINGES 2 PIECES INJEKT



## MAIN SPECIFICATIONS

- Composition: Medical grade plastic.
- Sterilized by ethylene oxide in packaging
- individual type "peel pack".
- Great quality.
- Clear and transparent cylinder.
- Clear range for accurate dosing.
- Safety washer.
- Smooth and regular feed of the plunger.
- Easy positioning of the plunger.
- Anatomical design of the plunger.

Code	Description	Selling unit
07001001	SERINGA 2 PEÇAS 2 ML LUER	CX.100
07001002	SERINGA 2 PEÇAS 5 ML LUER	CX.100
07001003	SERINGA 2 PEÇAS 10 ML LUER	CX.100
07001004	SERINGA 2 PEÇAS 20 ML LUER	CX.100



[www.factormed.pt](http://www.factormed.pt)